

Translation of the original operating instructions

Ultrasonic Med M, XL

Ultrasonic cleaning unit



RUCK[®] Ultrasonic Med M 2603602 RUCK[®] Ultrasonic Med XL 2603603 Copyright © 2023 Elma Schmidbauer GmbH. All rights reserved.

Manufacturer: Elma Schmidbauer GmbH Gottlieb-Daimler-Strasse 17 78224 Singen (Germany) www.elma-ultrasonic.com

Service and complaints: Tel. +49 7082 944 22 55 Email: service@hellmut-ruck.de Distribution: Hellmut Ruck GmbH Daimlerstrasse 23 75305 Neuenbürg Germany www.hellmut-ruck.de

Service Switzerland: Elma Schmidbauer Suisse AG Hintermättlistrasse 11 5506 Mägenwil Switzerland info@elma-suisse.ch www.elma-suisse.ch

Subject to technical and visual modifications.

Contents

1	Abou	ut this manual	4
	1.1	User	4
	1.2	Other applicable documents	4
	1.3	Factory settings	4
	1.4	CE mark	5
2	Prop	er use	5
	2.1	Cleaning medical devices	5
	2.2	Misuse	6
3	Safe	ty	6
	3.1	The structure of warnings	7
	3.2	Meaning of symbols on the unit	7
	3.3	Meaning of symbols on the nameplate	_
	. .		7
	3.4	Safety information on specific types of hazard	8
4	Use	of cleaning agents 1	0
	4.1	Permissible cleaning agents 1	0
	4.2	Cleaning agents 1	0
5	Tech	nical data 1	1
6	Prod	luct contents 1	3
7	Unit	description1	4
8	Cont	rols 1	5
	8.1	Ultrasound mode display 1	6
	8.2		6
	8.3		17
9	Set u	ip and connect the unit 1	17
•	9.1	-	17
	9.2		8
	- · · -		-

	9.3	Connect the unit	18
10	Fill th	ne tank	19
	10.1	Dosing cleaning agent	20
11	Inser	t objects to be cleaned	20
12	Oper	ating the unit	21
	12.1	Switch on the unit	21
	12.2	Set the cleaning duration	22
	12.3	Set the target temperature	22
	12.4	Set the ultrasound mode	22
	12.5	Set auto-degas	23
	12.6	Set Temperature-controlled start	23
	12.7	Start ultrasonic cleaning	23
	12.8	Save programmes	24
	12.9	Start a saved ultrasonic cleaning programme	24
	12.10	Switch off and save	25
	12.11	Remove cleaned objects	25
13	Drain	the unit	26
14	Chan	ge parameters in Setup	27
15	Malfu	inctions/fault messages	28
	15.1	Troubleshooting	28
	15.2	Fault messages	29
16	Main	tenance	30
	16.1	Replace the network cable	30
	16.2	Clean enclosure	30
	16.3	Cleaning the tank	30
	16.4	Disinfect	31
	16.5	LED test	31
17	Dispo	osal	31



1 About this manual

IMPORTANT Please read the manual before using the product. This manual forms part of the contents supplied with the unit. Keep the manual in an accessible place close to the product, and keep it with the product if it is sold on.

Meaning of the symbols used:

- This symbol denotes lists.
- This symbol denotes requirements.
- 1. Numbers with a dot denote actions.
- This symbol denotes individual actions.
- This symbol denotes interim results.
- → This symbol denotes the result of an action.
- 1 Numbers without a dot denote image labels.

1.1 User

In the manual, the term User refers to all persons who transport, set up, connect, operate and maintain the unit. The manual is aimed at persons with appropriate technical knowledge and experience in handling similar units.

The user is at least 16 years old. They must have read and understood the manual and be capable of following all notes and instructions.

All tasks that go beyond the pure operation of the unit within the scope described here must be performed by qualified and authorised specialists.

1.2 Other applicable documents

In addition to the specifications given in this manual, the following documents and directives also apply and may take precedence. This list is not exhaustive:

- Regional safety regulations and other locally applicable provisions.
- Safety data sheet and dosing instructions of the cleaning agents used.
- All serious incidents that occur in connection with the product must be reported to the manufacturer and the responsible authority of the member state in which the operator and/or the patient are resident.

1.3 Factory settings

The unit is delivered with predefined factory settings that determine the type and version of various unit functions. All information given in this manual is based on the use of a unit that is being operated with the factory settings unchanged.

Unit functions can be customised in *Setup* mode. If this has been done, the unit functions can differ from the information given in this manual. If discrepancies are found, the unit can be reset to the factory settings [▶ 27].

1.4 CE mark

This unit fulfils CE marking requirements in line with EU (EC/EEC) directives. Details are stated in the EU declaration of conformity that can be obtained from the manufacturer.

Technical changes to the unit must be approved by the manufacturer. The unit is no longer CE compliant if opened by an unauthorised person.

2 Proper use

The unit is used exclusively to treat objects and materials submerged in cleaning fluid with ultrasound. The cleaning effect is produced by the purely physical action of the ultrasound generated by the unit. The cleaning power can be boosted by adding suitable cleaning agents [▶ 10] and by increasing the temperature and duration of the bath. The unit is also authorised for use in typical laboratory applications, such as the preparation of samples by degassing, emulsifying or dispersing fluids and solutions.

The unit is classified as a class I medical device in line with the EU Medical Device Regulation and is designed exclusively for use in commercial environments. Only original accessories manufactured for the unit are permitted to be used.

The following objects are considered to be appropriate cleaning goods as long as the object in question is approved for cleaning in an ultrasonic bath. The list is not exhaustive:

- Medical and surgical microinstruments/instruments.
- Other medical devices.
- Dental instruments and dental laboratory materials.
- Podiatry instruments.
- Instruments from the tattooing and piercing sector.
- · Laboratory tools such as flasks, filters and plastic objects.
- Spectacles and objects used in timepiece and jewellery production.
- Industrial products and industrial objects.

Using the unit in any other way is considered improper use.

The operator is responsible for assessing the cleaning results.

2.1 Cleaning medical devices

The unit is used to clean and pre-clean medical devices. The following conditions must be fulfilled:

- The medical device is approved for ultrasonic cleaning and reprocessing (see information from the medical device manufacturer in line with EN ISO 17664).
- The medical device is only reused if downstream cleaning and preparation procedures ensure that it is properly cleaned, disinfected and sterilised.
- In cases of doubt, information on using and cleaning the device that is stipulated in the associated manual must be given precedence.

IMPORTANT Pre-cleaning medical products does not replace subsequent cleaning, disinfection or sterilisation using automated standard processes, e.g. washer-disinfectors or autoclaves.



3 Safety

2.2 Misuse

Misuse refers to any use of the unit that deviates from proper use. Misuse has the following consequences:

- Any misuse at own risk.
- The warranty is invalidated if the unit is misused.
- All liability for personal injuries and damage to property is excluded.

Therefore, it is vital to prevent misuse. This particularly applies to foreseeable misuse that can, for instance, arise in the following circumstances:

- Operation by persons who are incapable of operating the unit safely due to their physical, sensory or mental abilities, due to inexperience or due to a lack of knowledge. Children, in particular, must be kept away from the unit at all times.
- Failure to observe safety and warning notifications and maintenance and repair regulations.
- Failure to observe the conditions and fluids defined in the manual during setup and operation.
- Operation using unapproved cleaning agents.
- · Operation in areas with explosive atmospheres.

• Operation connected to a power supply without a fault-current circuit breaker.

Use for the following purposes is also considered misuse:

- · Cleaning or hydrating contact lenses.
- Cleaning animals or plants.

3 Safety

The unit has been designed and manufactured in line with current technological standards and recognised technical safety regulations. Nevertheless, the unit can pose a hazard to the life and health of users or third parties throughout its entire service life and can cause damage to itself or other material assets.

This chapter provides information on the basic structure of warnings, possible residual risks when using the unit and on how to avoid the resulting hazards. This information must be supplemented by statutory provisions and locally applicable regulations, which can be very different depending on the place of installation and intended application.

3.1 The structure of warnings



SIGNALWORD

Type and source of the danger

Possible consequences of the danger if disregarded.

> Measures to prevent the danger.

Signal word	Meaning	Consequences if disregarded		
DANGER	Immediate danger	Death or serious injury		
WARNING	Possible dangerous situation	Death or serious injury		
CAUTION	Possible dangerous situation	Slight injuries, damage to components or units		
NOTICE	Useful advice or tip	No risk of personal injury, but possible damage to components or units		

Table 1: Meaning of the signal words

3.2 Meaning of symbols on the unit

|--|

3.3 Meaning of symbols on the nameplate

\triangle	Attention
MD	Medical device
Ĩ	Read the manual
	Manufacturer
$\sim \sim$	Date of production
REF	Order number
SN	Serial number
-15°C	Temperature limit
CE	CE mark



3 Safety





Disposal information

3.4 Safety information on specific types of hazard

Electrical power

Contact with live components can lead to serious

injury or death due to electrocution!

- Do not connect the unit to the power supply if the connection cables or other components are visibly damaged.
- · Keep the enclosure and control elements clean and dry.
- Protect the unit against penetrating moisture.
- The mains voltage and the connected load on the nameplate must correspond to the on-site connection conditions.
- Only operate the unit in power circuits that are protected by a fault-current circuit breaker.

Risk of fire and explosion

Critical injuries, burns

- Do not fill the cleaning tank with flammable liquids.
- · Only use cleaning agents that are approved for this unit.
- · In case of doubt, ask the manufacturer or supplier.

Infections

Danger of possible infection due to poor cleaning results, dirty cleaning fluid and insufficient maintenance and disinfection of the unit.

- Change the cleaning fluid when it becomes visibly dirty and at least once per day.
- Clean the tank and surfaces thoroughly after draining the fluid and disinfect, if necessary.
- Check the cleaning results. The operator is responsible for monitoring the cleaning results.

Hot liquids and surfaces

Danger of burns and scalds due to contact with hot fluids or surfaces caused by high operating temperatures or continuous ultrasound operation.

- Do not touch the surfaces, accessories or objects being cleaned.
- Danger of splashes caused by high temperatures, by switching on ultrasound operation or by careless insertion of the basket or objects. If necessary, wear suitable protective equipment.
- If work must be performed on hot components, switch off the unit and allow to cool. If necessary, wear suitable protective equipment.

Cleaning agents

Volatile, corrosive or aggressive cleaning agents can cause chemical burns to the skin and respiratory tract.

- Observe the safety data sheet when using cleaning agents.
- Wear the protective equipment stipulated in the safety data sheet.
- Where necessary, provide adequate extraction for the emitted vapours and regularly check that the extraction system is working correctly.
- Observe the information stipulated in the *Cleaning agents* chapter. In case of doubt, ask the manufacturer or supplier.

Ultrasound-conductive liquids and materials

Ultrasound damages the cell membranes and bone structure.

- Do not reach into the ultrasonic bath during ultrasound operation.
- During ultrasound operation, do not touch any parts that conduct ultrasound, such as the tank, basket or any accessories used.

Ultrasound noise emissions

Working with ultrasonic cleaning units for longer periods of time can damage your hearing.

- Use the lid or wear personal ear protection when working on ultrasonic units.
- Pregnant women must not be subjected to the noise emissions for longer periods of time.
- · Keep animals away from the vicinity of ultrasonic devices.

Electromagnetic radiation

• No binding statement can be made for persons with active medical implants, such as pacemakers or implanted defibrillators. A binding statement can only be made for specific workstations and in consultation with the implant manufacturer.



4 Use of cleaning agents

4 Use of cleaning agents

A cleaning agent can be added to the liquid in the tank to improve the cleaning results.

If using cleaning agents, observe and apply the instructions stated in the safety data sheet and product information.

Also observe the following information:

4.1 Permissible cleaning agents

- Only use aqueous cleaning agents for ultrasonic cleaning.
- Only use cleaning agents that are suitable for ultrasonic cleaning.

4.2 Cleaning agents

Use of aggressive, corrosive cleaning agents:

Aggressive, corrosive cleaning agents can cause severe chemical burns if they come into contact with eyes or skin. Insufficient ventilation can result in severe chemical burns to the respiratory tract.

- Before using any cleaning agents, read the safety data sheet and product information and observe and apply the instructions.
- Where necessary, install an effective extractor system for the emitted vapours (e.g. at higher temperatures) and regularly check that the extractor system is working correctly.
- · Never reach into the tank when filled with cleaning agent.
- Wear suitable protective equipment in line with the safety data sheet.

Danger of fire or explosion when using combustible cleaning agents

Ultrasound and heat increase the evaporation of liquids and form extremely fine mists that readily ignite on contact with sources of ignition. This can result in severe burns or death.

- Do not use cleaning agents that are labelled with the pictograms GHS01 (explosive), GHS02 (flammable) or GHS03 (oxidising) in line with the CLP regulation (EC no. 1272/2008). Do not use cleaning agents that have a flashpoint.
- If necessary, clarify the cleaning agents that can be used by consulting the manufacturer or supplier.



Table 2: GHS pictograms for explosive, flammable or oxidising substances

Damage to stainless-steel tank caused by unsuitable cleaning agents

Unsuitable cleaning agents can cause pitting and therefore damage the stainless-steel tank within a short time.

- Consult the unit and cleaning agent manufacturers to clarify whether the cleaning agent is suitable.
- Do not use any cleaning agents in the acidic pH range at the same time as halides, such as fluoride, chloride, bromide or iodide.
- Only use cleaning agents that are suitable for ultrasound applications.

5 Technical data

Ultrasonic Med	Unit	М	XL	
Mechanical data				
Max. external dimensions W/D/H (unit with cover)	mm	305 / 170 / 230	365 / 275 / 280	
Max. opening dimensions W/D (inside of tank above stacking collar)	mm	235 / 135	295 / 235	
Interior dimensions W/D (tank below fill level marking)	mm	220 / 120	280 / 220	
Fill level (tank base to fill level marking)	mm	65	115	
Working height (basket base to fill level marking)	mm	40	80	
Interior dimensions of basket W/D/H (original accessories)	mm	195 / 105 / 50	255 / 200 / 70	
Max. basket loading (original accessories)	kg	1	6	
Total volume	I	2,7	9,5	
Recommended tank operating capacity	I	1,6	6,8	
Weight (unit with cover)	kg	3,6	6,2	
Drain connection thread	Inches	-	3/8	
Nominal hose nozzle diameter	-	-	DN 10	
Material (tank, enclosure)	-	Stainless	steel	
Performance data				
Total power consumption	W	320	550	
Ultrasonic power effective ±10 %	W	120	150	
Ultrasonic peak performance max.	W	480	600	
Heating power	W	200	400	
Electrical data				
Mains voltage ±10%	V~	220 - 240		
Mains frequency	Hz	50 / 6	60	
Ultrasonic frequency -2.5 / +5.5	kHz	37		
Max. target temperature 25 – 80	°C	30 (Temperature tolerance max. from -5 - +		
Degree of protection	_	IP 20		
Protection class	-	I		



5 Technical data

Ambient conditions					
Temperature (transport)	°C	-15 - +60			
Temperature (operation, storage)	°C	+5 - +40			
Air pressure (transport, storage)	hPa	500 - 1010			
Permissible relative humidity (transport,	% r. h.	10 – 80, non-condensing			
Permissible relative humidity (operation)	% r. h.	80; non-condensing under fluctuating temperatures			
Max. perm. altitude (operation)	m above sea level	+2000			
Overvoltage category	-	II			
Degree of contamination	-	2			
Sound pressure level LpAU*	dB	< 80			

*Measured sound pressure level with cover at a distance of 1 m.

6 Product contents

Check deliveries for damage to the packaging. Document any damage immediately

(e.g. photo), and report it to the manufacturer or dealer.

IMPORTANTCheck that all parts of the delivery are complete and undamaged.
Never put a damaged unit into operation.

Dispose of packaging materials that are no longer required in an environmentally friendly manner.



Illustration 1: Product contents

1	Ultrasonic unit Ultrasonic Med M / XL
2	Cover or, turned over, usable as a drip tray for a basket
3	Basket
4	Hose connector for drain (Ultrasonic Med XL, not pictured)
5	Network cable (IEC connector, not shown in the image)
6	Manual (not shown)



7 Unit description

7 Unit description



Illustration 2: Operating side/rear side

- 1 Display, splash-proof
- 2 Warning of hot surfaces
- 3 Handle
- 4 Nameplate
- 5 Mains connection
- 6 Outflow (Ultrasonic Med XL)
- 7 Discharge tap, outflow open/closed (Ultrasonic Med XL)
- 8 Membrane keyboard, wipeable
- 9 Feet, non-slip

8 Controls




Illustration 3: Controls

- 1 Cleaning duration display value/unit/symbol
- 2 Current ultrasound mode display name/bar segment
- 3 **Start/stop button** for ultrasonic cleaning (short press)/ultrasonic cleaning with auto-degas (long press)
- 4 Temperature display value/unit/symbol
- 5 **Plus/Minus buttons** set temperature/program memory/parameter values (short press); quick cycle (long press)
- 6 Heating button on/off (short press)/temperature-controlled start (long press)
- 7 **On/off button** for switching the unit on/off
- 8 **Mode button** set ultrasound modes (short press)/switch modes (long press)/save programs (press > 4 seconds)
- 9 Plus/minus buttons set cleaning duration/program memory/parameters
- 10 **Display** showing all LED indicators for LED test LED test [> 31]



8 Controls

8.1 Ultrasound mode display

The display shows the currently set ultrasound mode by displaying the mode name and a certain number of bar segments. You can choose from the following ultrasound modes:

Biga.	degas (0 segments) auto-degas (degas flashing)	Used to degas cleaning fluid and for specialist laboratory applications. Used to degas the cleaning fluid before ultrasonic cleaning. The set ultrasound mode starts automatically once degassing has been completed (duration approx. 10 min.).
	eco (1 segment)	Used for gentle cleaning and quieter operation. The <i>sweep</i> ultrasound mode is a lower-intensity cycle.
SUCCS	<i>sweep</i> (2 segments)	This mode performs most cleaning tasks with an evenly distributed sound field (standard ultrasound mode).
pulse	<i>pulse</i> (3 segments)	Removes stubborn dirt using pulsed sound field distribution.
dyn.	<i>dynamic</i> (4 segments)	Combines the ultrasound modes <i>sweep</i> and <i>pulse</i> to improve the overall cleaning power.

8.2 Operating status symbols

The following symbols are shown on the display as soon as the unit enters the respective operating status:

~	End of cleaning	Ultrasonic cleaning has been completed and is finished.
٢	Remaining time	Ultrasonic cleaning is running; the remaining cleaning duration is shown below the symbol.
	Dry operation	(Not available for Elmasonic Med 500/900). The tank is empty, or the filling level is too low. Dry-run detection is available during active heating.
\$	Setup	Setup mode is activated; the unit parameters [▶ 27] can be changed.
ß	Heating	The ultrasonic bath is currently being heated.
J o	Temperature- controlled start	Ultrasonic cleaning will start automatically as soon as the set cleaning temperature has been reached.
≬ ▼	Temperature limit	A temperature limit is defined as a maximum permissible cleaning temperature. A higher temperature cannot be set. The symbol flashes when the temperature limit is reached or exceeded.
1	Programme	The displayed cleaning programme is selected (programme number 1–4). A total of four cleaning programmes can be stored.

8.3 Operating modes

Operating mode	Brief description	Display		
Switched off	The unit is being supplied with power and is switched off.	Off		
Standby	The unit is switched on. The display is dimmed after a longer period of inactivity. The automatic safety shut-off is activated after 8 hours of inactivity.	On		
Setup	The unit is switched off. The Setup menu is activated to change the unit functions.	 Setup symbol Number of the selected unit function (parameter) Current setting (parameter value) 		
Normal mode	The unit is ready to use. Ultrasonic cleaning can be conducted using adjustable operating parameters.	Set operating parametersDuring Normal mode:Current operating status		
Programme mode	The unit is ready to use. Four stored programmes can be selected for ultrasonic cleaning. The programmes can be specified and saved.	 Available programmes Operating parameters of the selected programme No programmes have been saved: Programme number 1 and "" for duration and temperature During programme operation: Current operating status 		

9 Set up and connect the unit

9.1 Notes on installation site

The unit is designed to be set up in commercially used premises such as laboratories, medical practices, etc.

- The installation site must be well ventilated.
- The permissible ambient conditions are met.
- The unit must be installed on a sturdy, horizontal, non-slip, moisture-resistant work surface with the following characteristics:
 - The surface must be non-flammable.
 - There must be sufficient clearance above the unit to allow the cover and the objects being cleaned to be removed.
 - There must be an area available next to the unit for setting down the basket and inserts.
- Electrical connection conditions:
 - Splash-proof earthed socket near the unit (cable length approx. 1.5 m).
 - The power connection must be protected by a fault-current circuit breaker (RCD).
 - Ensure that the power supply required to operate the unit is available.



9 Set up and connect the unit

- Wastewater connection requirements:
 - Units without drains: select an installation site located close to a sink that can be used to empty the unit.
 - Units with drains: properly connect the unit drain to the on-site wastewater disposal system. If necessary, properly route the drain into a sink located nearby.
 CAUTION! Use heat-resistant material.

9.2 Set up the unit

- ✓ Unpack the unit completely.
- ✓ Keep the cover and accessories e.g., basket, inserts, etc., on hand.
- Place the unit on the selected surface.
 IMPORTANT The controls, such as the display, buttons and, if applicable, discharge tap must be clearly visible and easily accessible.
- → The unit has been set up.

9.3 Connect the unit

Connect the drain

- ✓ Obtain the required connection materials, e.g. suitable heat-resistant hose and hose clips (the hose is not included).
- 1. Unscrew the yellow plug from the drain.
- 2. Screw the included hose nozzle into the drain.
- 3. Push the hose onto the hose nozzle and secure it using the hose clip.
- 4. Connect the hose to the wastewater disposal system. If necessary, properly route the hose into a sink.
- 5. Check the drain for leaks. To do so, fill the tank with cold water up to around 1/3 of the fill level. Then, open the discharge tap to test all joints for leaks.
- 6. Seal any leaks found.
- ➔ The drain is connected.

Connecting the network cable

- ✓ Ensure that the required power connection is available.
- 1. Plug the power cord into the IEC connector on the unit.
- 2. Route the power cord in such a way that it does not pose a tripping hazard, cannot be damaged and is not exposed to moisture. The power cord must not touch the unit housing, as the housing becomes hot during operation.
- 3. Connect the plug. The plug must be easily accessible so that it can be disconnected in emergencies.
- → The power supply has been connected.

10 Fill the tank

Suitable cleaning liquids for ultrasonic cleaning are:

- Water
- Softened water
- Distilled water

Cleaning power can be improved by adding diluted cleaning agents. See Dosing cleaning agent [> 20].

- ✓ The unit and the tank are clean and, if necessary, disinfected.
- ✓ The recommended tank operating capacity for the relevant unit sizes is stated in the technical data. See Technical data.
- ✓ The unit is ready to use.
- 1. Add the fluid to the tank up to the fill level marking using, for example, a measuring jug. The fill level marking is the stacking collar (1). For unit sizes without stacking collars, the fill level is approx. 1.5 cm beneath the tank edge.

CAUTION! An overfilled tank can cause the cleaning fluid to boil over at high temperatures.



- 2. The process of adding cleaning agents is explained in chapter Dosing cleaning agent [> 20].
- Fresh liquid must be degassed approx. 10 minutes prior to ultrasonic cleaning. See Setting degas: Set the ultrasound mode [▶ 22] or Set auto-degas [▶ 23].
 IMPORTANT! It is important to degas freshly inserted ultrasonic cleaning baths to remove any air bubbles and dissolved gases that may be present in the cleaning liquid. Failure to remove these gases will negatively impact cleaning performance.
 - \Rightarrow The tank has been filled.

Insert a basket:

1. Hang the basket filled with the items to be cleaned in the tank.

IMPORTANT! The fill level in the tank will increase once the basket or insert is placed inside. This may cause the tank to overflow in the case of full baskets or dosages > 5%. To avoid this, do not fill the tank right up to the fill level marking in such cases and fill the remaining space once the basket has already been inserted.

- 2. **IMPORTANT!** The objects being cleaned must be completely submerged in the fluid to achieve good cleaning results. If necessary, reduce the number of objects being cleaned.
- ➔ The tank is full, cleaning agent has been dosed (if necessary), liquid has been degassed and the basket has been hung. Ultrasonic cleaning can be started.



11 Insert objects to be cleaned

10.1 Dosing cleaning agent

Cleaning power can be improved by adding diluted cleaning agents.

A WARNING! Please follow the safety instructions and instructions for use and observe the mixing ratio for the cleaning agent used!

Dosing information for the different unit sizes and typical volume concentrations are shown in the following table.

Unit size	Tank operating capacity	Dosing quantities in mI for following concentrations:			owing
Volume concentration in %		1%	2%	3%	5%
	1 l/1000 ml	10 ml	20 ml	30 ml	50 ml
Ultrasonic Med S	0.7 l / 700 ml	7 ml	14 ml	21 ml	35 ml
Ultrasonic Med M	1.6 l / 1600 ml	16 ml	32 ml	48 ml	80 ml
Ultrasonic Med XL	6.8 l / 6800 ml	68 ml	136 ml	204 ml	340 ml

✓ The diluted cleaning agent is chosen based on the type of cleaning desired.
 ▲ WARNING! Do not use flammable cleaning agents! See Use of cleaning agents [▶ 10].

- ✓ The unit is filled with liquid (e.g. water) up to the recommended tank operating capacity. See Fill the tank [▶ 19].
- 1. Find the dosing quantity for the required concentration in the table and add the liquid to the unit.
- 2. Stir the cleaning liquid to mix. Stirring is also achieved by degassing the cleaning liquid. See Setting degas: Set the ultrasound mode [▶ 22] or Set auto-degas [▶ 23].
- → The cleaning agent is dosed into the liquid. For further steps, see Fill the tank [▶ 19].

11 Insert objects to be cleaned

- Ensure that all optional original accessories e.g., baskets or inserts for holding the objects being cleaned, are available.
 IMPORTANT Do not place objects on the bottom of the tank. The tank and the objects can be damaged through ultrasonic operation.
- 1. Position the basket or insert next to the unit. The cover can be used as a mat.
- Place the objects inside the basket or insert without overfilling it. Arrange the objects being cleaned in such a way that fluid will be able to circulate around them without obstructions.
 IMPORTANT Only clean objects that are suitable for ultrasonic operation, the set temperature and, where applicable, the cleaning agent in use. In case of doubt, contact the manufacturer or dealer.
- 3. Hook the filled basket or insert inside the tank.
- → The objects being cleaned have now been inserted.

12 Operating the unit

WARNING



Hearing damage!

- When working with the unit in ultrasound mode, cover it with the lid or wear hearing protection.
- Pregnant women must not spend long periods of time near an ultrasonic unit in operation.
- > Keep animals away from the ultrasonic unit.

Damage caused by noise emissions

CAUTION



Burns or scalds.

Do not reach into the tank.

Hot liquid, vapours and unit parts!

When working at high cleaning temperatures (> 50 °C), only touch baskets, covers and inserts wearing suitable protective gloves or allow them to cool before touching them.

CAUTION



Dry running!

Damage to the tank and unit.

- > Always fill the tank with fluid up to the fill level marking before operating the unit.
- > Never operate the unit without cleaning fluid.
- > Regularly monitor the filling level.

12.1 Switch on the unit

The unit is ready to use.

The unit is filled and the objects to be cleaned have been inserted.

The unit is covered with the cover.




Press the **On/Off** button to switch on the unit. The display shows the last-used settings, e.g. duration, ultrasound mode, current temperature and, if applicable, the stored programme number.




12 Operating the unit


Two modes can be selected for starting ultrasound operation:

Display in Normal mode



- 1. Set the duration, ultrasound mode and temperature.
- 2. Switch the heating on/off.
- 3. Set the unit to start with auto-degas.
- 4. Or, set the unit to start with temperature control.





Select or start a saved programme.

5. Save programmes.

Press the *Mode* button (1) for 2 seconds to switch between Programme mode and Normal mode. The display (2) shows the current settings.

The display (3) shows the set programme number (programmes 1–4 in Programme mode).

12.2 Set the cleaning duration

From 1 minute to 6 hours.



Press one of the *Plus/Minus* buttons below the time display to set the duration. Press the button repeatedly until the required duration is shown on the display. Hold the button down to increase/decrease the duration more quickly.

12.3 Set the target temperature

From 25-80 °C/77-176 °F.

Press one of the *Plus/minus* buttons below the temperature display to set the required target temperature. Press the button repeatedly until the required temperature flashes on the display. Hold the button down to increase/decrease the duration more quickly.

+

IMPORTANT If a temperature limit between 40–60 °C/104–140 °F has been set, see the Change parameters section for information on changing the temperature limit. Longer periods of ultrasound operation heat the cleaning fluid in addition to the heating function. This can lead to the fluid being hotter than the set target temperature.

12.4 Set the ultrasound mode

Degas, eco, sweep, pulse or dyn (dynamic)



Press the *Mode* button to set the ultrasound mode. Press the button repeatedly until the required ultrasound mode is shown on the display.

12.5 Set auto-degas

Use this function to degas a fresh batch of cleaning fluid.

One of the ultrasound modes, eco, sweep, pulse or dyn, has been set.

The unit is not currently in ultrasound operation.

Hold down the *Start/Stop* button until the ultrasound mode *degas* flashes on the display.

Press the Start/Stop button again to start the auto-degas function.

Alternatively, save auto-degas as a programme. See the section Save programmes for more information.

Degas flashes slowly on the display (degassing takes approx. 10 min.). Once auto-degas has been completed, the set ultrasonic cleaning mode will start and be shown on the display.

12.6 Set Temperature-controlled start

Use this function to start ultrasound operation automatically once the target temperature has been reached.

IMPORTANT When cleaning sensitive objects, please note that the ultrasound is switched on briefly multiple times during the heating phase.

The required target temperature has been set.

Ensure that the required target temperature is higher than the current actual temperature. If not, allow the fluid to cool.

The unit is not currently in ultrasound operation.



Hold down the *Heating* button until the *Temperature-controlled start* symbol flashes on the display.



Press the *Start/Stop* button to start the unit using the temperature-controlled start function.

See the section Save programmes for information on saving the Temperature-controlled start function as a programme.

The display shows the Temperature-controlled start and Heating symbols. The display shows the set target temperature and the current actual temperature of the cleaning fluid flashing alternately.



During the heating phase, the ultrasound is switched on briefly multiple times to ensure that the fluid is heated evenly.

The ultrasonic cleaning function will start once the set target temperature has been reached. The Temperature-controlled start symbol will disappear.

12.7 Start ultrasonic cleaning

The display shows the current settings.

For information on changing the settings, see Set cleaning duration, ultrasound mode and target temperature.

✓ The described settings have been configured and should be started.



If necessary, press the *Heating* button to heat the cleaning fluid. The Heating symbol will be shown on the display. The display shows the set target temperature and the current temperature of the cleaning fluid flashing alternately. If necessary, wait until the target temperature has been reached.



12 Operating the unit

Press the *Start/stop* button to begin ultrasound operation in the displayed mode. The Remaining time symbol is shown on the display. The remaining time, the current temperature and the set ultrasound mode are continuously shown on the display.



The ultrasound mode, temperature, duration and heating **on/off** can be changed during ultrasonic cleaning. **Important** Changing the duration switches off ultrasonic cleaning. Press the **Start/stop** button to restart ultrasonic cleaning with the altered time setting.




The End of cleaning symbol and, where applicable, an acoustic signal (depending on the configured parameters) indicate that cleaning has ended. The Remaining time display will disappear.

12.8 Save programmes

You can save up to 4 programmes.

Display in Save mode



Display in Programme mode



To save the currently set programme, use the plus/minus buttons to select a programme number from 1 to 4.

✓ In Normal mode the settings were made in the preset order. Switch on the unit [▶ 21] These settings can be saved as a programme.

Important The auto-degas and temperature-controlled start modes can only be saved separately as individual programmes.



Hold down the *Mode* button (1) for >4 seconds until the display shows all 4 programme numbers (3) and programme number 1 flashes. Now, press one of the *Plus/Minus* buttons repeatedly until the programme number under which you wish to save your settings flashes. Press the *Mode* button (1) again briefly to confirm saving the programme. The saved programme number (4) will be shown on the display.

12.9 Start a saved ultrasonic cleaning programme

✓ At least one programme from a possible four programme numbers has been saved.

The unit is in Programme mode. The last-used programme number is shown on the display.

✓ The unit is in Normal mode. No programme number is shown on the display.

If necessary, switch between Normal mode and Programme mode. See Switch on the unit [> 21].



Hold down the *Mode* button until a programme number is shown on the display. The last-used programme number is always shown, e.g. programme number 1.



Press one of the *Plus/minus* buttons to select a programme.



Press the *Start/stop* button to start the displayed programme.

The Remaining time and Heating symbols will be shown on the display (if saved). The remaining time, the current temperature and the ultrasound mode are continuously shown on the display.



IMPORTANT Depending on the settings saved in the programme, ultrasonic cleaning is started with auto-degas (degas flashes) or Temperature-controlled start (only starts once the target temperature has been reached). The respective symbols are shown on the display.



The end of cleaning is indicated by the *End of cleaning* symbol and, where applicable, an acoustic signal (depending on the parameter settings). The Remaining time display will disappear.

12.10 Switch off and save

For information on saving programmes, see Save programmes [> 24].

✓ The displayed settings are saved when the unit is switched off and are then available the next time the unit is switched on. This function cannot be used to save *Auto-degas* and *Temperature-controlled start*.



Press the **On/off** button to switch off the unit. Any settings made will be saved until the next time the unit is switched on.

12.11 Remove cleaned objects

- ✓ The end of ultrasonic cleaning is indicated by the End of cleaning symbol or, where applicable, an acoustic signal.
- 1. Carefully remove the lid, and allow any dripping water to drip into the tank. Then, place the lid upside down next to the unit to use as a drip tray.

A CAUTION! If using high cleaning temperatures > 50 °C, allow the liquid to cool down first or use suitable heat-insulating protective gloves.

2. Lift the basket or insert out of the tank, allow the excess fluid to drain briefly, and place it on the upturned lid.

⇒The objects are now clean. Check the cleaning result.

- → Check the fill level if you intend to start another ultrasonic cleaning process.
- ➔ Switch off the unit. If the fluid is dirty or will not be used again, empty out the liquid. See the section Draining the unit for information.



13 Drain the unit



CAUTION

Hot liquids

Danger of burns and scalds caused by spilled fluid!

- \succ Do not move or carry the unit when full.
- Exception: Units without drains. Allow the fluid to cool beforehand, and disconnect the mains plug.

Empty units with drains:

- ✓ Ultrasound operation has finished.
- ✓ The basket or insert has been removed.
- ✓ The drain (1) is connected to an on-site wastewater disposal system.
- Open the discharge tap (2) to drain off the cleaning fluid.
 ⇒The tank is empty.
- 2. Clean and, if necessary, disinfect the tank. See chapter Maintenance for more information. **CAUTION! Never submerge the unit in water.**
- → The unit is emptied, cleaned and, where necessary, disinfected.



Empty units without drains:

- ✓ Ultrasound operation has finished.
- ✓ The basket or insert has been removed.
- ✓ If using temperatures > 50 °C, allow the cleaning fluid to cool.
- 1. Disconnect the network cable.
- 2. Put the cover on.
- 3. Firmly take hold of the unit by the handles while holding the cover in place, and carry it carefully to a sink.
- 4. Remove the cover.
- 5. Slowly tip the unit over the sink with one corner at the lowest point, and carefully pour out the cleaning fluid.

⇒The tank is empty.

- 6. Clean and, if necessary, disinfect the tank. See chapter Maintenance for more information. **CAUTION! Never submerge the unit in water.**
- → The unit is emptied, cleaned and, where necessary, disinfected.

14 Change parameters in Setup

IMPORTANT If parameters marked with an asterisk (*) are changed, the saved programme will be deleted. If necessary, make a note of the settings of the saved programmes before changing them.

- \checkmark The unit is connected to a power supply.
- ✓ The Parameter reference list is required for setting the parameter numbers and values.
- The unit has been switched off.



Hold down the *Mode* button until the Setup menu is shown on the display.



First, press the *Plus/minus* buttons to set the parameter number, e.g., 01.



Then, press the *Plus/minus* buttons to change the parameter value, e.g., 00/01. Set and change further parameters as required.



Press the **On/off** button to save the settings.

Parameter reference:

The following unit functions can be adjusted. The preset factory values are shown in bold:

Parameter	Unit function	Value	Effect
01 *	Adjustment range for cleaning duration > 10 min.	00 /01	Increments of 5 min./increments of 1 min. (as for cleaning durations up to 10 min.). NOTICE! Saved programmes will be deleted.
02	Show remaining time < 1 min. in seconds	00/ 01	Off/ on
03 *	Temperature unit	°C /°F	Display in degrees Celsius / Fahrenheit. NOTICE! Saved programmes will be deleted.
04 *	Default temperature limit	00 /40/45/50/55/60 00 /100/110/120/ 130/140	Off/40 °C/45 °C/50 °C/55 °C/ 60 °C Off/100 °F/110 °F/120 °F/130 °F/ 140 °F NOTICE! Saved programmes will be deleted.
05	Acoustic signal for "End of cleaning"	00/ 01 /02	Off/melody 1/melody 2
06	Acoustic signal for "Error message"	00/ 01 /02	Off/melody 1/melody 2

15 Malfunctions/fault messages

07	Acoustic signal for "Temperature limit exceeded"	00/ 01 /02	Off/melody 1/melody 2
08	Response if temperature limit is exceeded	00/ 01 /02/03	No action/ ultrasound off /heating off/ultrasound and heating off
09	Display brightness	01/02/ 03 /04	Dim/medium/normal/bright
10	LED test	00 /01	Off /all display LEDs on for approx. 5 sec.
20 ***	Autostart on restoration of power supply	00 /11/12/13	Off/ultrasound on/heating on/ ultrasound and heating on CAUTION! Risk of damage to the unit if the unit runs dry.
30 *	Reset unit settings	00 /11	Off/on: "11" resets all settings to the factory settings. NOTICE! Saved programmes will be deleted.
31 **	Response if dry-run monitor is triggered	00/ 01	No action/ ultrasound and heating off, display Dry run symbol.

** Not available for unit sizes 500/900.

*** The unit starts up immediately in heating and/or ultrasound operation when plugged in or when the power supply is switched on.

15 Malfunctions/fault messages

15.1 Troubleshooting

Fault	Possible cause	Resolution
Display does not light up.	Unit not switched on or no power available.	 Switch on the unit. Check that the power cable is correctly connected and undamaged; replace if necessary. Check the power supply.
Buttons not working.	Controls defective.	Contact the manufacturer.
Not possible to set a high bath temperature.	A low temperature limit has been defined.	 Check the unit settings [> 27]; increase or deactivate the temperature limit as required.
Displayed bath temperature incorrect.	Temperature sensor defective.	Contact the manufacturer.
Heating not working.	Heating element defective.	Contact the manufacturer.

Fault	Possible cause	Resolution
Ultrasonic operation does not start.	Ultrasonic unit defective.	Contact the manufacturer.
Display shows <i>Er</i> .	Error message.	Check the error code and rectify the fault; contact the manufacturer if necessary.
Blower is off (only units with blowers)	Ultrasound is off and bath temperature < 50 °C or blower defective.	 Check the blower if the ultrasound is switched on or the bath temperature is > 50 °C; contact the manufacturer if necessary.

15.2 Fault messages

If a fault occurs, the display shows the type of fault using a symbol or error number. You can rectify the following faults yourself:

Display	Fault type	Possible remedies
	Filling level too low (dry run detected).	Top up fluid. If the symbol continues to be displayed, have the heating system repaired.
) •	Temperature limit exceeded.	Switch off the unit and allow it to cool. Disconnect from mains power if necessary.
Er 21	Max. perm. bath temperature exceeded.	Switch off the unit and allow it to cool. Disconnect from mains power if necessary.
Er 40–44	Power supply faulty.	Connect the unit to a socket of a different power circuit.
Er 45 -46	Ultrasonic power faulty.	Change fill level, reduce/increase if needed Remove basket, if necessary, reduce the number of objects being cleaned.

Table 3: Error messages

Contact the service centre or the manufacturer in the following cases:

- The measures described here do not remedy the faults.
- The faults recur (repeatedly), even after performing the following measures:
 - Disconnect the unit from the mains power, then reconnect and restart the unit after approx.
 1 minute.
 - The unit has been reset to factory settings.
- The error number displayed is not listed here: Make a note of the error number; the unit must be repaired at the factory.

IMPORTANT Empty the unit completely and clean and disinfect it thoroughly before sending it for repair.

16 Maintenance



Faulty power cord

CAUTION

Electric shock or damage to the device

> Check the power cord regularly for damage.

IMPORTANT Keep the unit and tank clean and dry to prolong its service life. Regularly remove all dirt residues and accumulated dust from the blower.

16.1 Replace the network cable

- ✓ The network cable or plug contacts are visibly damaged.
- Replace damaged network cables immediately.
 NOTICE! Removable network cables may not be replaced with insufficiently measured network cables!
- → The network cable has been replaced.

16.2 Clean enclosure

- \checkmark The enclosure is dirty.
- ✓ The network cable is unplugged.
- Wipe all surfaces using a damp cloth.
 WARNING! Do not spray the outside of the unit with water.
- → The enclosure has been cleaned.

16.3 Cleaning the tank

- ✓ The tank is full of limescale or dirty.
- ✓ The tank is empty.
- 1. Use a damp cloth to wipe out limescale and dirt residue. If necessary, rinse the tank using a shower set.

CAUTION! Do not spray the outside of the unit with water.

⇒The tank has been cleaned.

- 2. To descale the tank, fill it with tap water up to the fill level marking in the tank.
- 3. Heat up the unit to 40 °C.
- 4. Add 4–10% of Elma clean 60 or Elma clean 115C.
- 5. Once 40 °C has been reached, switch off the unit.

 \Rightarrow Leave the solution to work for approx. 12 hours.

- 6. Then, switch on eco ultrasound mode for approx. 15 min.
- 7. Empty the tank.
- 8. Use a damp cloth to wipe out the remaining limescale and dirt residue. If necessary, rinse the tank using a shower set.

CAUTION! Do not spray the outside of the unit with water.

→ The tank has been descaled and cleaned.

16.4 Disinfect

We recommend the following disinfectants for wiping the unit:

- 29666 RUCK ® alcohol spray and wipe disinfection.
- The unit is used in the medical and health sector.
 NOTICE! Observe local directives and cleaning regulations. This particularly applies to the cleaning of medical units.
- Disinfect the tank and the housing using a commercially available surface disinfectant on a regular basis. A CAUTION! Test the disinfectant on a small area first to ensure it does not affect the materials, particularly the controls.
- → The unit has been hygienically disinfected.

16.5 LED test

Test the LED function (see Controls)

Regularly check the function of all LEDs on the display to ensure reliable operation. The display showing all LED indicators is illustrated in the chapter Controls. Send the unit for repair immediately if any LED indicators are not displayed.

The unit is switched off.



Hold down the *Mode* button until the Setup menu is shown on the display.



To set the LED test function, press the *Plus/minus* buttons repeatedly until the parameter number 10 is shown.



Then, press the *Plus* button to change the parameter value from 00 to 01.



The LEDs will light up for approx. 5 seconds. The parameter value will be reset automatically. Repeat this procedure as required.

 $\bigcirc \bigcirc \bigcirc$

Press the **On/off** button to end the LED test.

17 Disposal

CAUTION

Once the unit has reached the end of its service life, ensure that the unit and accessories are disposed of safely and correctly:



- Clean and disinfect the old device and accessories before disposal.
- Do not dispose of old devices with household waste, but instead at the local collection and disposal points.
- Secure the old device against unauthorised access until removal; if necessary, dispose of the power cable separately.
- > Observe regionally applicable disposal directives.
- Data protection notice: The end user is responsible for deleting personal and confidential data from the unit being discarded.



HELLMUT RUCK GmbH Daimlerstraße 23 DE–75305 Neuenbürg fon +49 (0)7082 944 20 fax +49 (0)7082 944 22 22 kontakt@hellmut-ruck.de hellmut-ruck.de