

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 1 of 11

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

Baehr Desinfektionsmittelkonzentrat für Instrumente

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Disinfectant for processing of instruments
Product for professional use**1.3. Details of the supplier of the safety data sheet**

Company name:	Gustav Baehr GmbH	
Street:	Max-Eyth-Str. 39	
Place:	71332 Waiblingen	
Telephone:	+49 (0) 7151 95902-0	Telefax: +49 (0) 7151 95902-62
e-mail:	info@gustav-baehr.de	
Responsible Department:	Responsible for the safety data sheet: sds@gbk-ingelheim.de	

1.4. Emergency telephone number:INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a)
In England and Wales: NHS 111 In Scotland: NHS 24 - dial 111**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Regulation (EC) No. 1272/2008**Hazard categories:
Skin corrosion/irritation: Skin Irrit. 2
Serious eye damage/eye irritation: Eye Dam. 1
Hazardous to the aquatic environment: Aquatic Acute 1
Hazardous to the aquatic environment: Aquatic Chronic 2
Hazard Statements:
Causes skin irritation.
Causes serious eye damage.
Very toxic to aquatic life.
Toxic to aquatic life with long lasting effects.**2.2. Label elements****Regulation (EC) No. 1272/2008****Hazard components for labelling**N,N-didecyl-N-methyl-poly(oxyethyl)ammoniumpropionate
Alcohol ethoxylate C9-C11
Guanidine, N,N"-1,3-propanediylbis-, N-coco alkyl derivs., diacetates
Tridecylamine**Signal word:** Danger**Pictograms:****Hazard statements**

H315	Causes skin irritation.
H318	Causes serious eye damage.
H410	Very toxic to aquatic life with long lasting effects.

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 2 of 11

Precautionary statements

P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P311	IF exposed or concerned: Call a POISON CENTER/doctor.
P501	Dispose of contents/container to approved disposal company or local collection.

2.3. Other hazards

According to Regulation (EC) No 1907/2006 (REACH) none of the substances, contained in this product are a PBT / vPvB substance.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterization

Aqueous preparation with detergents and solvents

Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	GHS Classification			
112-34-5	2-(2-Butoxyethoxy)ethanol			< 5 %
	203-961-6	603-096-00-8	01-2119475104-44	
	Eye Irrit. 2; H319			
169107-21-5	Alcohols, C9-11-branched, ethoxylated			< 5 %
	Acute Tox. 4, Eye Dam. 1; H302 H318			
94667-33-1	N,N-didecyl-N-methyl-poly(oxyethyl)ammoniumpropionate			< 5 %
	619-057-3		01-2119950327-36	
	Acute Tox. 4, Skin Corr. 1B, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1; H302 H314 H318 H400 H410			
85681-60-3	Guanidine, N,N"-1,3-propanediylbis-, N-coco alkyl derivs., diacetates			< 5 %
	288-198-7		01-2119980967-14	
	Acute Tox. 4, Skin Corr. 1C, Eye Dam. 1, Aquatic Acute 1, Aquatic Chronic 1; H302 H314 H318 H400 H410			
106232-83-1	Alcohols, C12-15-branched and linear, ethoxylated			< 5 %
	Acute Tox. 4, Eye Dam. 1, Aquatic Chronic 3; H302 H318 H412			
107-21-1	Ethane-1,2-diol			< 2,5 %
	203-473-3	603-027-00-1	01-2119456816-28	
	Acute Tox. 4, STOT RE 2; H302 H373			
86089-17-0	Tridecylamine			< 2,5 %
	289-185-9		01-2119461722-40	
	Acute Tox. 4, Skin Corr. 1B, STOT SE 3, Aquatic Acute 1, Aquatic Chronic 1; H302 H314 H335 H400 H410			

Full text of H and EUH statements: see section 16.

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 3 of 11

Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
		Specific Conc. Limits, M-factors and ATE	
112-34-5	203-961-6	2-(2-Butoxyethoxy)ethanol	< 5 %
		dermal: LD50 = 4120 mg/kg; oral: LD50 = 5660 mg/kg	
169107-21-5		Alcohols, C9-11-branched, ethoxylated	< 5 %
		oral: ATE = 500 mg/kg	
94667-33-1	619-057-3	N,N-didecyl-N-methyl-poly(oxyethyl)ammoniumpropionate	< 5 %
		oral: LD50 = 1157 mg/kg M akut; H400: M=10 M chron.; H410: M=1	
85681-60-3	288-198-7	Guanidine, N,N"-1,3-propanediylbis-, N-coco alkyl derivs., diacetates	< 5 %
		oral: ATE = 500 mg/kg M akut; H400: M=10 M chron.; H410: M=1	
106232-83-1		Alcohols, C12-15-branched and linear, ethoxylated	< 5 %
		oral: ATE = 500 mg/kg	
107-21-1	203-473-3	Ethane-1,2-diol	< 2,5 %
		dermal: LD50 = 9530 mg/kg; oral: LD50 = 5840 mg/kg	
86089-17-0	289-185-9	Tridecylamine	< 2,5 %
		oral: LD50 = 820 mg/kg M akut; H400: M=10 M chron.; H410: M=10	

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

Remove contaminated soaked clothing immediately.
If you feel unwell, seek medical advice.

After inhalation

Move to fresh air in case of accidental inhalation of vapours or decomposition products.
In the event of symptoms refer for medical treatment.

After contact with skin

Wash off immediately with soap and plenty of water.
Consult a doctor if skin irritation persists.

After contact with eyes

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes.
Consult (eye) doctor immediately.

After ingestion

Do not provoke vomiting. Consult physician. Attention in case of vomiting - acute danger of suffocating, produced by foaming ingredients. Rinse mouth. Make drink some glasses of water. The decision whether to provoke vomiting is to be taken by a physician.

4.2. Most important symptoms and effects, both acute and delayed

Causes skin irritation.
Causes serious eye damage.
Watch out. Beware, hazard of foam aspiration.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptoms.

SECTION 5: Firefighting measures

5.1. Extinguishing media

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 4 of 11

Suitable extinguishing media

Product does not burn, fire-extinguishing activities according to surrounding.

Unsuitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixtureFire may produce:
carbon monoxide (CO), carbon dioxide (CO₂) and nitrogen oxides (NO_x)**5.3. Advice for firefighters**

Wear self-contained breathing apparatus and protective suit.

Additional information

Cool containers at risk with water spray jet.

Fire residues and contaminated firefighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General measures**Ensure adequate ventilation.
Use personal protective clothing.**For non-emergency personnel**Do not breathe vapours.
Avoid contact with skin, eyes and clothing.**For emergency responders**In case of vapour formation use respirator.
Use personal protective clothing.**6.2. Environmental precautions**

Do not discharge into the drains/surface waters/ground water.

6.3. Methods and material for containment and cleaning up**For containment**

Prevent spread over a wide area (e.g. by containment or oil barriers).

For cleaning upSoak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder).
Shovel into suitable container for disposal.**6.4. Reference to other sections**Observe protective instructions (see Sections 7 and 8).
Information for disposal see section 13.**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**Keep container tightly closed.
Avoid contact with skin, eyes and clothing.
Use only in well-ventilated areas.**Advice on protection against fire and explosion**

No special protective measures against fire required.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container tightly closed in a dry, cool and well-ventilated place.

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 5 of 11

Hints on joint storage

Incompatible with oxidizing agents.

Further information on storage conditions

Keep away from food, drink and animal feeding stuffs.

7.3. Specific end use(s)

Disinfectant for processing of instruments

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Exposure limits (EH40)**

CAS No	Substance	ppm	mg/m ³	fibres/ml	Category	Origin
112-34-5	2-(2-Butoxyethoxy)ethanol	10	67.5		TWA (8 h)	WEL
		15	101.2		STEL (15 min)	WEL
107-21-1	Ethane-1,2-diol, vapour	20	52		TWA (8 h)	WEL
		40	104		STEL (15 min)	WEL

8.2. Exposure controls**Appropriate engineering controls**

Ensure adequate ventilation, especially in confined areas.

Protective and hygiene measures

Wash hands before breaks and immediately after handling the product.

When using do not eat, drink or smoke.

Take off immediately all contaminated clothing.

Avoid contact with skin, eyes and clothing.

Eye/face protection

Eye wash bottle with pure water (EN 15154).

Tightly fitting goggles (EN 166).

Safety goggles with side protection (EN 166).

Hand protection

Chemical-resistant gloves (EN 374)

(recommended: minimum protection index 2, corresponding to a permeation rate > 30 minutes according to EN 374)

Butyl rubber - Butyl (0,5 mm)

Fluoro-rubber - FKM (0.4 mm)

This recommendation refers exclusively to the chemical compatibility and the lab test conforming to EN 374 carried out under lab conditions.

Requirements can vary as a function of the use. Therefore it is necessary to adhere additionally to the recommendations given by the manufacturer of protective gloves.

Skin protection

Long sleeved clothing (DIN EN ISO 6530)

Respiratory protection

In case of insufficient ventilation wear suitable respiratory equipment (gas filter type A) (EN 14387).

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state:

Liquid

Colour:

Yellowish

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 6 of 11

Odour: Perfumed (scented)

Test method

pH-Value: 8,0 - 10,0 Concentrate

Changes in the physical state

Melting point: n.d.

Boiling point or initial boiling point and boiling range: approx. 100 °C

Sublimation point: n.a.

Softening point: n.d.

Flash point: n.a.

Sustaining combustion: Not sustaining combustion

Flammability

Solid/liquid: n.a.

Gas: n.a.

Explosive properties

The product is not explosive.

Lower explosion limits: n.d.

Upper explosion limits: n.d.

Auto-ignition temperature: n.d.

Self-ignition temperature

Solid: n.a.

Gas: n.a.

Decomposition temperature: n.d.

Oxidizing properties

Not oxidising.

Vapour pressure: n.d.

Density (at 20 °C): 0,992 - 1,012 g/cm³

Bulk density: n.a.

Water solubility: Miscible
(at 20 °C)**Solubility in other solvents**

n.d.

Partition coefficient n-octanol/water: n.d.

Viscosity / dynamic:
(at 20 °C) n.d.

Viscosity / kinematic: n.d.

Flow time: n.d.

Relative vapour density: n.d.

Evaporation rate: n.d.

Solvent separation test: n.d.

Solvent content: < 10 %

9.2. Other information

No data available

SECTION 10: Stability and reactivity

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 7 of 11

10.1. Reactivity

No decomposition if stored and applied as directed.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Reactions with oxidizing agents.

10.4. Conditions to avoid

To avoid thermal decomposition, do not overheat.

10.5. Incompatible materials

Strong oxidizing agents

10.6. Hazardous decomposition products

No hazardous decomposition products known.

Fire may produce:

Carbon monoxide (CO), carbon dioxide (CO₂) and nitrogen oxides (NO_x)**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute toxicity**

Based on available data, the classification criteria are not met.

No toxicological data available.

ATEmix/oral: > 2000 mg/kg

Irritation and corrosivity

Causes skin irritation.

Causes serious eye damage.

Skin: Non corrosive [OECD 435]

Sensitising effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

Classification in compliance with the assessment procedure specified in the Regulation (EC) no 1272/2008 and on the basis of toxicological data.

11.2. Information on other hazards**Other information**

Ingestion causes irritation of upper respiratory system and gastrointestinal disturbance.

Watch out. Beware, hazard of foam aspiration.

SECTION 12: Ecological information**12.1. Toxicity**

Ecological data are not available.

Very toxic to aquatic life with long lasting effects.

12.2. Persistence and degradability

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 8 of 11

No data available

12.3. Bioaccumulative potential

No data available

12.4. Mobility in soil

No data available

12.5. Results of PBT and vPvB assessment

According to Regulation (EC) No 1907/2006 (REACH) none of the substances, contained in this product are a PBT / vPvB substance.

12.7. Other adverse effects

Severe hazard to waters

Further information

Do not flush into surface water or sanitary sewer system.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**Where possible recycling is preferred to disposal.
Can be incinerated, when in compliance with local regulations.**List of Wastes Code - residues/unused products**

070699 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of fats, grease, soaps, detergents, disinfectants and cosmetics; wastes not otherwise specified

Contaminated packagingEmpty containers should be taken for local recycling, recovery or waste disposal.
Contaminated packaging should be emptied as far as possible.
Packaging that cannot be cleaned should be disposed of like the product.**SECTION 14: Transport information****Land transport (ADR/RID)**

14.1. UN number:	UN 3082
14.2. UN proper shipping name:	ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. (Disinfectant)
14.3. Transport hazard class(es):	9
14.4. Packing group:	III
Hazard label:	9



Classification code:	M6
Special Provisions:	274 335 375 601
Limited quantity:	5 L / 30 kg
Excepted quantity:	E1
Transport category:	3
Hazard No:	90
Tunnel restriction code:	-

Other applicable information (land transport)

Special provision 375: Product in packaging (inner packaging, single packaging) up to 5 L is not subject to the provisions of ADR/RID.

Inland waterways transport (ADN)

according to Regulation (EC) No 1907/2006


Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 9 of 11

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
(Disinfectant)
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9



Classification code: M6
Special Provisions: 274 335 375 601
Limited quantity: 5 L / 30 kg
Excepted quantity: E1

Other applicable information (inland waterways transport)

Special provision 375: Product in packaging (inner packaging, single packaging) up to 5 L is not subject to the provisions of ADN.

Marine transport (IMDG)

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
(Disinfectant)
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9



Marine pollutant: Yes
Special Provisions: 274, 335, 969
Limited quantity: 5 L / 30 kg
Excepted quantity: E1
EmS: F-A, S-F

Other applicable information (marine transport)

Subsection 2.10.2.7 IMDG Code: Product in packaging (inner packaging, single packaging) up to 5 L is not subject to the provisions of MDG Code.

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
(Disinfectant)
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9



Special Provisions: A97 A158 A197
Limited quantity Passenger: 30 kg G
Passenger LQ: Y964
Excepted quantity: E1

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 10 of 11

IATA-packing instructions - Passenger:	964
IATA-max. quantity - Passenger:	450 L
IATA-packing instructions - Cargo:	964
IATA-max. quantity - Cargo:	450 L

Other applicable information (air transport)

Special provision A197: Product in packaging (inner packaging, single packaging) up to 5 L is not subject to the provisions of IATA/DGR.

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: Yes

**14.6. Special precautions for user**

Take the usual precautions when handling with chemicals.

14.7. Maritime transport in bulk according to IMO instruments

The transport takes place only in approved and appropriate packaging.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3, Entry 55

2004/42/EC (VOC): 7 %

Information according to 2012/18/EU (SEVESO III): E1 Hazardous to the Aquatic Environment

National regulatory information

Employment restrictions:

Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.

Water hazard class (D):

3 - highly hazardous to water

15.2. Chemical safety assessment

For this substance a chemical safety assessment has not been carried out.

SECTION 16: Other information**Changes**

Changes in section: 2, 3, 11, 12, 14, 15

Abbreviations and acronyms

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route

RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses

ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure

IMDG = International Maritime Code for Dangerous Goods

IATA/ICAO = International Air Transport Association / International Civil Aviation Organization

MARPOL = International Convention for the Prevention of Pollution from Ships

IBC-Code = International Code for the Construction and Equipment of Ships Carrying Dangerous Chemicals in Bulk

according to Regulation (EC) No 1907/2006

Baehr Desinfektionsmittelkonzentrat für Instrumente

Revision date: 30.04.2021

Product code: 00320-0027N-GHS

Page 11 of 11

GHS = Globally Harmonized System of Classification and Labelling of Chemicals
REACH = Registration, Evaluation, Authorization and Restriction of Chemicals
CAS = Chemical Abstract Service
EN = European norm
ISO = International Organization for Standardization
DIN = Deutsche Industrie Norm
PBT = Persistent Bioaccumulative and Toxic
vPvB = Very Persistent and very Bio-accumulative

LD = Lethal dose
LC = Lethal concentration
EC = Effect concentration
IC = Median immobilisation concentration or median inhibitory concentration

Classification for mixtures and used evaluation method according to Regulation (EC) No. 1272/2008 [CLP]

Classification	Classification procedure
Skin Irrit. 2; H315	Bridging principle "Dilution"
Eye Dam. 1; H318	Calculation method
Aquatic Acute 1; H400	Calculation method
Aquatic Chronic 2; H411	Calculation method

Relevant H and EUH statements (number and full text)

H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.
H315 Causes skin irritation.
H318 Causes serious eye damage.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.
H373 May cause damage to organs through prolonged or repeated exposure.
H400 Very toxic to aquatic life.
H410 Very toxic to aquatic life with long lasting effects.
H411 Toxic to aquatic life with long lasting effects.
H412 Harmful to aquatic life with long lasting effects.

Further Information

Data of items 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities.

The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge.

The delivery specifications are contained in the corresponding product sheet.

This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

(n.a. = not applicable; n.d. = not determined)

'The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)